**ergolet**

# Walking slings

*Walking sling – male  
Walking sling – female*

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## Important

This manual must be read before the sling is used. This manual is intended for therapists and other personnel, in connection with the use and maintenance of *Ergo* slings.

## Purpose

Ergolet's *Ergo* slings can be used for transfers in home care, nursing homes, institutions and private homes. *Ergo* slings are used to transfer people with physical or reduced mobility issues.

## Check list before use

It is the carer's responsibility to check the sling before use.

- It must be ensured that the sling does not have any visible damages or wear and tear and that the sling is clean.
- When the sling is used, qualified personnel must be present.
- The capacity written on the sling must never be exceeded!

## Lifetime

The sling can be washed 80-100 times, which corresponds to once a month in 6-8 years. The product's lifetime is expected to be 5-10 years, if service checks and the stated maintenance instructions have been followed. It is the owner's responsibility that the annual inspection is performed.

## Materials

The slings are produced in durable materials and tolerate the necessary washing. All slings are produced in polyester.

## Label symbols

Declaration of Conformity

Recommended 92°C for washing

Do not bleach

Do not iron

Do not dry clean

Tumble dry at low temperature

Max 180kg

Maximum lifting capacity is 180kg

## Placement of sling straps

### *Colour coded straps*

The straps on most of the slings are marked with colours (red, yellow, green, blue), which makes adjustment of position and placement of the sling easier.

Corresponding sling straps must be attached to the spreader bar with the same colour. The sitting position can be varied by choosing different colours for the leg and shoulder straps.

### *Patented Ergolet sling hook system*


Ergolet does not take any responsibility if the sling is placed incorrectly. It is always the carer's responsibility that the correct sling has been chosen and the sling is fitted correctly.

1. The sling strap is placed on the open hook (1).
2. The strap is pulled downwards and the hook will close (2).

The design of the lifting hook is such that the weight of the user will help close the hook.


To release the strap after lifting simply pull upwards to open the hook.

Note: Ensure that client's fingers cannot get trapped.

### *Strap color code*

RED	SMALL
YELLOW	MEDIUM
GREEN	LARGE
BLUE	X-LARGE

## Lifting

Walking slings for gait training with either Ergolet Luna Ceiling Hoist or Weight Relief System.

Note: The sling may NEVER be used to lift the user. The sling is only for safety purposes.

A qualified instructor or carer must always be present during gait training with the Weight Relief System.

The sling can be attached to the user while sitting in a chair (Fig. 1 and Fig. 2).

1. Place the sling in the chair and tighten the belts around the user, while making sure that the sling fits.
2. Attach leg strap with buckle and place the spreader bar straps on the spreader bar.
3. Lift the user up until standing on both legs.

NB: The user must never be lifted by this sling while standing.

- Ensure that sling loops are placed correctly on the spreader bar hooks.
- Ensure the patient sits/lies safely before taking off the sling.

Above: A user in the Walking sling (male version) attached on the Luna ceiling hoist.

Above: Gait training on a treadmill in the Walking sling (female version).

# Walking slings

## USABILITY

For gait training e.g. rehabilitation. When the user needs extra safety.

Available in female and male versions.

Walking sling  
Female model



Walking sling  
Male model

## TECHNICAL SPECIFICATION

Max capacity	180 kg
Max washing temperature	92° C
Max drying temperature	55° C
Spreader bar compatibility	2-point
Sling is available in sizes	XXS - XL

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